



EVERFLOW SCIENTIFIC INSTRUMENTS

Old #298, New #231, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.
Ph: 26447066, 26447708, Fax: 91-044-26447708 Grams: EVERFLOW
web site: everflowscientificinstruments.com, e-mail: everflow@eth.net

VERTICAL AUTOCLAVE



SPECIFICATIONS:

- Everflow vertical autoflow is electrically heated steam sterilizer is used for sterilization of surgical instruments, dressing material, lien, rubber, etc.
- Working Pressure 15 to 20 psi.
- Sturdy double walled construction with inner chamber made of Stainless Steel. Outer Chamber made of Mild Steel Shet duly powder coated.
- Lid is made of Stainless Steel Plate and is tightened alround by wing nuts.
- Moulded jointless gaskets are made of Neoprene rubber.
- ISI marked Immersion heating elements heats the water and steam to desired temperature and pressure.
- The Autoclave is provided with Pressure Gauge, Steam release Valve, Spring loaded safety Valve.
- Work on 220/230 Volts, A.C.

SIZES:

	H	x	D
a)	450	x	250 mm
b)	500	x	300 mm
c)	550	x	350 mm
d)	600	x	400 mm
e)	600	x	450 mm
f)	750	x	550 mm
g)	650	x	450mm.

OPTIONS AT EXTRA COST

- a) Automatic Low Water Cut Off Device.
- b) Pressure Control Switch.
- c) Digital Temperature Controller cum Indicator with timer.

DRY STEAM & HIGH PRESSURE AUTOCLAVES ARE ALSO AVAILABLE

Mfrs. : ★ Scientific ★ Surgical ★ Hospital ★ Laboratory ★ R & D ★ Hydrology ★ Meteriologic
★ Survey ★ Industrial Instruments ★ Educational Instruments ★ Equipments