


BALL MILL



For mixing & grinding, electrically operated having capacity of few gms to 1 KG & 2 KG. Fitted with geared motor with max speed of 80 rpm. The jar is made of SS and having steel balls of different sizes to be put in jar for mixing & grinding pesticides power. Jar is interchangeable and easily removable, light in weight and compact unit working on 220 V AC with 1 KG cap.

- do -- SS Drum 2 KG capacity
- do -- Aluminium Drum 1 Kg Capacity
- do -- Aluminium drum 2 KG capacity