


DELUXE MODEL

Hand Operated Centrifuge
12/24 Butyrometers
Heavy Duty Centrifuge

SUPER DELUXE CENTRIFUGE HAND OPERATED:

FOR 24 BUTYROMETERS
HEAVY DUTY CENTRIFUGE
WITH PROTECTIVE COVER



Accessories for
use in 'Gerber'
fat test



Electric Gerber Centrifuge

Dual Purpose For 8/12 Butyrometers
Or 4/6 Solubility Index Tube.
With Timer 0-5 Minutes, Mechanical
Brake & Anti Vibration Rubber Pads

ELECTRIC GERBER CENTRIFUGE

24 Tests. Midi Super Model.
Protective Cover, With Timer
0-5 Minutes, Mechanical Brake.
Anti Vibration Rubber Pads



EVERFLOW SCIENTIFIC INSTRUMENTS

No. 231/298, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.
M : 9444259624, LL : 044-26646066, e-mail : everflowin@gmail.com