

MEMBRANE FILTRATION ASSEMBLY



Membrane Filter Holder to hold 47mm. dia. Filter

- Stainless Steel Membrane filter holder
(OR)
- Glass Membrane Filter holder

Filter Flasks: Capacity : 1000ml.

OIL-FREE VACUUM PUMPS: Diaphragm type

Model No.	Max. Flow (Ltrs/Min.)	Max. Vacuum.	Max. Pressure (Psig)	Motor HP
OFD-15-S	15	22"	25	1/20
OFD-15-D	15	27"	35	1/16

Filter paper Suitable for the above

Cellulose Nitrate Filter Paper, 47mm dia, 0.45 microns.
Cellulose Nitrate Filter Paper, 47mm. dia. 0.2 microns.

www.everflowscientific.com



EVERFLOW SCIENTIFIC INSTRUMENTS

No. 231 / 298, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.

Mob: 9444259624, Ph : 044-26646066, 26645708 e-mail: everflowin@gmail.com

GSTIN: 33AFBPG5986P1ZX