



EVERFLOW SCIENTIFIC INSTRUMENTS

Old #298, New #231, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.

Ph: 26646066, 26645708, Tele Fax: 91-044-26645708, Mobile No: 94442 59624

E-Mail: everflow@eth.net, Web Site: www.everflowscientific.com

HOT AIR OVEN



SPECIFICATIONS:

- ▶ Everflow Oven is Double walled construction with glass wool insulation.
- ▶ This Outer chamber made of Mild Steel sheet duly powder coated.
- ▶ Inner Chamber made of Stainless Steel Sheet.
- ▶ Temperature Range from ambient to 250°C with an accuracy of $\pm 1^\circ\text{C}$ controlled by Thermostat.
- ▶ Heating elements made of high grade imported nichrome wire.
- ▶ Heating elements are placed in the Bottom for Lab Type and in three sides for Memmert Type.
- ▶ Door fitted with heavy casted chromium plated hinges with spring and automatic roller to lock.
- ▶ Works on 220/230 Volts A.C.

SIZES:

Available (Lab. Type/ Universal Type)	W	X	D	X	H	No. Of Shelves
a)	300	x	300	x	300	2
b)	355	x	355	x	355	2
c)	455	x	455	x	455	3
d)	455	x	455	x	605	3
e)	605	x	605	x	605	3
f)	605	x	455	x	910	3
g)	605	x	605	x	910	3

OPTIONS AT EXTRA COST

- a) Air Circulating Fan to obtain uniform temperature throughout the chamber.
- b) Electronic Digital Temperature Controller Cum Indicator instead of Thermostatic Control.
- c) Ovens upto 300 C Range are also available.
- d) Timer

OTHER PRODUCTS: VACCUM OVEN & INDUSTRIAL OVENS ETC.

Mfrs. : ★ Scientific ★ Surgical ★ Hospital ★ Laboratory ★ R & D ★ Hydrology ★ Meteorologic
 ★ Survey ★ Industrial Instruments ★ Educational Instruments ★ Equipments