

## KJELDAHL DIGESTION / DISTILLATION ASSEMBLIES

### Kjeldahl Digestion Assembly



KJEDHAL DIGESTION Assembly : 3 test of 500ml. flask Mantle Heaters with Indicator and Regulator for 3 test of 500ml. capacity Flask. Glass parts comprises of Kjeldahl Flask (500ml.) and Fume Glass Duct Tube with clamp and Boss Head.

### Kjeldahl Distillation Assembly



Kjeldahl Distillation Assembly : 3 test of 500ml. flask. Mantle Heaters with Indicator and Regulator for 3 test of 500ml. capacity Flask. Glass Parts comprises of : Kjeldahl Flask , Multiple adapter, two necks, Splash head, vertical, Liebig condenser, Receiver Adapter, separating funnel, Condenser Clamp with Boss head and Rubber Tube.

### KJELDAHL DIGESTION AND DISTILLATION ASSEMBLY COMBINED MODEL



Kjeldahl Digestion and Distillation Assembly supplied complete with The Combined assembly is mounted on a Stand made of MS Sheet with powder coating in an attractive colour. Kjeldahl digestion & distillation combine unit open element. Mantle type is suitable for 500ml, Flask with pilot lamp. Provided with condenser rack.

Glass parts Kjeldahl Distillation & Digestion Apparatus combined: (Set of 3 Nos.).

Distillation : Kjeldahl Flask, Multiple adapter, two necks, Splash head, vertical, Liebig condenser, Receiver Adapter, separating funnel, Condenser Clamp with Boss head and Rubber Tube.

Digestion : Kjeldahl Flask, Condenser Clamp with Boss Head and Glass Fume Duct Tube .

EVERFLOW SCIENTIFIC INSTRUMENTS

No 231/298, Kilpauk Garden Road, Kilpauk, Chennai – 600 010.

Ph : 044-26646066, e-mail : [everflowin@gmail.com](mailto:everflowin@gmail.com) web : [www.everflowscientific.com](http://www.everflowscientific.com)