



EVERFLOW SCIENTIFIC INSTRUMENTS

Old #298, New #231, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.
 Ph: 26447066, 26447708, Fax: 91-044-26447708 Grams: EVERFLOW
 web site: everflowscientificinstruments.com, e-mail: everflow@eth.net

LAMINAR AIR-FLOW



VERTICAL LAMINAR AIR-FLOW

SPECIFICATIONS:

- Everflow Laminar Air-flow are constructed by using best quality Industrial block boards. Exteriors covered by decorative laminated sheets. Inner side coated with synthetic rubber paint.
- HEPA FILTERS are constructed from non-woven fibre glass media and zero probed to achieve better than 99.99% efficiency on particulate matter 0.3 microns; with greater efficiency on larger or smaller particles.
- Exclusive design results in a...uniform velocity over the sterile work zone.
- Airflow system prevents uneven particle loading by eliminating direct blasts to the HEPA filter spreading a uniform static pressure over 100% of the filter surface.
- HEPA filter protective guard helps to eliminate damages to the HEPA filter by external sources.
- Powder coated, well balanced uniform speed, blower assembly ensures constant flow.
- Standard accessories - illumination, castor wheels, Manometer, power point and gas valve.

HORIZONTAL LAMINAR FLOW				
Working Area	Overall Dimension (mm)			
	FT	L	D	H
2' x 2' x 2'	680	1050	1460	
3' x 2' x 2'	980	1050	1460	
4' x 2' x 2'	1280	1050	1460	
2' x 2' x 2'	1880	1050	1460	

VERTICAL LAMINAR FLOW				
Working Area	Overall Dimension (mm)			
	FT	L	D	H
2' x 2' x 2'	680	750	1900	
3' x 2' x 2'	980	750	1900	
4' x 2' x 2'	1280	750	1900	
2' x 2' x 2'	1880	750	1900	



HORIZONTAL LAMINAR AIR-FLOW

BIO SAFETY CABINET & FUME EXHAUST HOOD ARE ALSO AVAILABLE

Mfrs. : ★ Scientific ★ Surgical ★ Hospital ★ Laboratory ★ R & D ★ Hydrology ★ Meteriologic
 ★ Survey ★ Industrial Instruments ★ Educational Instruments ★ Equipments