

'EVERFLOW' MICRO SLIDE CABINET



'EVERFLOW' Micro Slide Cabinet

- The Micro slide Cabinet is used for storing Micro Slides of size 3"x1" in vertical position. Available in 5000 / 10000 slide capacities.
- Safe and Easy handling of slides.
- Cabinet made of teak ply wood body and the bottom trays of hard board with machine cut anodized aluminum groove, for holding slides.
- Slides are arranged separately in the tray without disturbances for other slides or compartments.
- The slide holding slotted rows are numbered to identify the slides.
- Complete with lacquer polish .
- The drawers will not be removed from the cabinet.
- The Cabinet is provided with a smooth operating door, with lock & key.

EVERFLOW SCIENTIFIC INSTRUMENTS

No. 298, Kilpauk Garden Road, Kilpauk, Chennai – 600010.

Phone No: 26646066, Tele fax : 044-26645708, website : www.everflowscientific.com

e-mail : everflowin@gmail.com, sales@everflowscientific.com