




M-TFW

M-TFW TISSUE FLOTATION WORKSTATION

It incorporates a feather touch keyboard with a LCD display. The temperature is controlled with the help of precise temperature controller. It also has a built in memory to store the set temperature and LED illumination inside the water bath for easy viewing. The water bath is made up of corrosion material which provides the required durability.

Technical Specifications

Water temperature range	: ambient to 100°C
Temperature Tolerance	: $\pm 1^\circ \text{C}$
Water bath dimensions, inner (mm)	: 258x193x64
Integrated hot plates	: Yes
Hot plate Dimensions (mm)	: 110x70 & 140x70
Integrated LED for bath illumination	: Yes



M-PD

M-PD - PARAFFIN DISPENSER:

It has a fully automatic procedure and memory function to store the set temperature, the feather touch keyboard and large display which displays the current heating state of the dispenser makes it easy to use.

The nozzle is fitted with temperature control so that the free flow of wax is maintained.

Technical Specifications

Range of Temperature	: ambient to 100°C
Capacity	: 5L/10L
Power	: 1200VA
Weight	: 17kg