

ULTRASONIC CLEANER BATH



Tank Capacity	Model No.	Tank Size (Moulded Tanks)	Ultrasonic power
1 Ltrs	MT - 1	150 x 135 x 65 mm (H)	50 watts
2 Ltrs	MT - 2	176 x 162 x 100 mm (H)	60 watts
3 Ltrs	MT - 3	240 x 135 x 100 mm (H)	80 watts
4.5 Ltrs	MT - 4	300 x 150 x 100 mm (H)	100 watts
6.5 Ltrs	MT - 6	300 x 150 x 150 mm (H)	150 watts
10 Ltrs	MT - 10	300 x 238 x 150 mm (H)	250 watts
14 Ltrs.	MT - 14	300 x 238 x 200 mm (H)	300 watts
20 Ltrs	MT- 20	500 x 300 x 150 mm (H)	400 watts
30 Ltrs.	MT - 30	500 x 300 x 200 mm (H)	500 watts

Timer	: 0-20 mins (For 1 & 2 ltrs) & 0 – 99 mins (for 3.3, 6.5,10,20 & 30)
Drain	: Will be available from 4.5 ltrs & onwards.
Basket	: Will be provided.
Heater	: temp cut off at 45° C (from 3.3 ltrs onwards)
U/s frequency	: 33 ± 3 KHz
Digital Temp Controller	: Will be provided

www.everflowscientific.com



EVERFLOW SCIENTIFIC INSTRUMENTS

No. 231 / 298, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.

Mob: 9444259624, Ph : 044-26646066, 26645708 e-mail: everflowin@gmail.com

GSTIN: 33AFBPG5986P1ZX