













New ATN-G SERIES

ANALYTICAL BALANCE - AUTO-CAL

AXIS

Manufactured in Poland



-  *animals weighing*
-  *density*
-  *force measurement*
-  *recipe*
-  *statistics*
-  *pieces counting*
-  *percentage*
-  *summation*
-  *weighing speed selection*
-  *threshold/selection*
-  *max weight*
-  *min weight*

Features

- * Complete automatic time & temperature calibration with GMP-GLP record printout.
- * Brilliant Graphic LCD display with simple navigation windows and menu to guide you through balance functions.
- * User selectable filter, stability indicator and auto-zero tracking for use in difficult environment.
- * Procell software allows for direct transfer of weighing results to excel spreadsheet.
- * Weigh below hook for below weighing application with in-built density measurement software.
- * With a die-cast metal base and stainless steel pan, and balance is built to withstand rugged every day Weighing.
- * Application : Parts Counting, Animal Weighing, Check-Weighing, Display Hold, Percent Weighing, Totalizing, SQC Quality Control.
- * RS232C and USB connector to connect with PC
- * Optional LAN / wi-fi

Option :-

-  *alibi memory*
-  *products and users database*
-  *printout config.*
-  *measurements database*

MODEL	CAPACITY (gm)	RESOLUTION (gm)	PAN (mm)	REPEATABILITY (gm)	LINEARITY (gm)	RESPONSE TIME (Sec.)
-------	---------------	-----------------	----------	--------------------	----------------	----------------------

Auto-Internal Calibration

ATN 120G	120	0.0001	∅ 90	± 0.0001	± 0.0002	< 3
ATN 220G	220	0.0001	∅ 90	± 0.0001	± 0.0002	< 3

Optional Accessories

- Density Kit • Data Printer • Software • Test Weights

