


‘EVERFLOW’ LAMINAR FLOW BENCH



Horizontal Laminar Flow Bench

Vertical Laminar Flow Bench

EVERFLOW’ LAMINAR FLOW CABINETS

EVERFLOW offers a comprehensive range of high performance Horizontal, Vertical & Bio Safe Laminar Flow Clean Air Cabinets. These are designed to meet the filtration, illumination, noise & vibration requirements of US Federal Standard providing particle free air to meet class 100 (class I) conditions.

- MOC : MS or Galvanized Iron Sheet with Powder Coating.
- Work Table : Made of SS 304 with side panels made of 6mm. thick transparent glass duly framed and Front door made of acrylic sheet.
- Filtration : Atmospheric air drawn through Pre-filter and pass through HEPA Filter having efficiency rating as high as 99.99% with cold DOP and 99.97% with hot DOP.
- Blower Assembly : Statically and dynamically balanced, direct drive, sized to provide adequate air flow volume over the entire surface of HEPA filter. Fitted with 1 Phase motor and operates with minimum noise level i.e., lower than 65 db on scale and Vibration less than 2.5 um. (0.0001 deg).
- Illumination : Work area properly illuminated by LED light.
- Types : Horizontal Laminar Flow Benches with Horizontal Flow of Filtered air.
: Vertical Laminar Flow Benches with Air Flows Downwards to Work Table.

Horizontal Laminar Flow Bench		Vertical Laminar Flow Bench	
<u>Model No.</u>	<u>Working Area (L x W x H)</u>	<u>Model No.</u>	<u>Working Area (L x W x H)</u>
EHL-22G	2’ x 2’ x 2’	EVL-22G	2’ x 2’ x 2’
EHL-32G	3’ x 2’ x 2’	EVL-32G	3’ x 2’ x 2’
EHL-42G	4’ x 2’ x 2’	EVL-42G	4’ x 2’ x 2’

Note : Laminar Flow Benches can also be constructed out of Commercial Ply wood or Complete SS based on user requirement at additional costs.

www.everflowsscientific.com



EVERFLOW SCIENTIFIC INSTRUMENTS

No. 231 / 298, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.

Mob: 9444259624, Ph : 044-26646066, 26645708 e-mail: everflowin@gmail.com

GSTIN: 33AFBPG5986P1ZX