

## ORBITAL SHAKER - Bench Top

A sturdy and reliable shaker, ideal for shaking flasks in all types of applications. Excellent for biotechnology.



- Body: MS painted with Epoxy powder, with SS tray.
- Drive: Permanent magnet geared DC motor.
- Shaking speed: 20-200 orbital strokes per min.
- Orbit dia.: 25mm.

### Available in following Platform capacities

Platform size : 325 x 325mm  
Clamps ; 100ml. – 16 Nos. (or) 250ml. – 8 Nos (or) 500ml. – 8 Nos.

Platform size : 425 x 425mm.  
Clamps suitable : 100ml. – 25Nos. (or) 250ml. – 16 Nos. (or) 500ml. – 16 Nos.

\*With Digital RPM Indicator.

\*Without Digital RPM Indicator.

[KINDLY SPECIFY THE CLAMPS WHILE ORDERING]

[www.everflowscientific.com](http://www.everflowscientific.com)



## **EVERFLOW SCIENTIFIC INSTRUMENTS**

No. 231 / 298, Kilpauk Garden Road, Kilpauk, Chennai - 600 010.

Mob: 9444259624, Ph : 044-26646066, 26645708 e-mail: [everflowin@gmail.com](mailto:everflowin@gmail.com)

GSTIN: 33AFBPG5986P1ZX