

## 'EVERFLOW' WATER BATHS

### WATER BATH RECTANGULAR – WITH CONCENTRIC RINGS (HOLES)



### 'EVERFLOW' WATER BATH – ( DOUBLE WALLED) : with MICROPROCESSOR BASED DIGITAL TEMPERATURE CONTROLLER CUM INDICATOR :

- Double walled Inner SS TIG Welded, Outer MS / SS.
- Glass wool insulated.
- Drain cock to facilitate easy emptying and cleaning of the bath.
- SS Lid with 3” Holes and covered by Concentric rings to hold 50 to 500ml. flasks.
- ISI marked U shaped Immersion heater to attain temperature of upto 99°C .
- Temperature controlled by Microprocessor based Digital Temperature Controller cum Indicator.

#### Available in following Inner Bath sizes

350 x 250 x 100 mm.	6 Holes	1000 watts
400 x 300 x 100 mm.	12 Holes	1500 watts



## **EVERFLOW SCIENTIFIC INSTRUMENTS**

Old No: 298, New No: 231, Kilpauk Garden Road, Kilpauk, Chennai - 10.  
Ph: 26646066, Fax: 26645708, Cell: 9444259624, Email: everflowin@gmail.com,  
TIN No: 33581122305, CST No: 866329 Dt: 27.06.2006