

'EVERFLOW' WATER DISTILLATION UNIT WALL MOUNTING TYPE



- Electrically operated.
- Made of Heavy gauge Stainless Steel of AISI 304 Grade.
- Lid is automatically sealed by water forming in a recess in the head of the Water Still into which the flange of lid rests.
- This minimizes the loss of steam.
- Still is fitted with Immersion type Heaters.
- Heating element is easily replaceable.
- Special care should be exercised to avoid the still running dry.

Supplied complete with Clamp for Wall Mount, chord & plug for operation on 230 volts 50 c/s. AC mains.

Available in following capacities

<u>Capacity</u>	<u>Rating</u>
2 Liters	2.0KW.
4 Liters	3.0 Kw
6 Liters	4.5Kw

EVERFLOW SCIENTIFIC INSTRUMENTS

No. 231/298, Kilpauk Garden Road, Kilpauk, Chennai-600010.
Chennai-600 010. Ph:044 26646066, Mob:9444259624
e-mail : everflowin@gmail.com, www.everflowscientific.com

